

Gamete Washing & Handling

PBS1X-Dulbecco's Phosphate Buffered Saline Solution

Cat #: APBS-500

500mL, PETG bottles

Description & Intended Use: Sterile Phosphate-Buffered saline solution (PBS1X) is intended for use in ART procedures such as washing gamete during retrieval/collection or labware washing prior to use. **Does not include Protein Supplement.**

Recommended Use: Ready-to-Use after equilibration at +37°C in incubator or ambient atmosphere based on the intended use. Consult with your laboratory protocols for additional specifics on the use of this product. Use aseptic techniques and discard any remaining media after the completion of the procedure. Use aseptic techniques.



- Lowest endotoxin levels in the market using most sensitive LAL testing <math><0.005\text{ EU/mL}</math>
- Triple QC extensive testing for each batch (Biocompatibility 1-Cell expanded MEA, LAL Endotoxin tested)
- Sterile Filtered and Sterility tested as per current USP <math><71></math> standards
- Manufactured under strict Quality Management System by ISO 13485:2016 certified facility
- FDA Registered Manufacturing Facility
- Batch-to-Batch consistency due to extensive quality control system and highly controlled process. COA for each released batch & SDS available upon request.
- Free sampling and free consultation on the use of the products is available upon request. Special pricing with standing orders and supply agreements.

Ordering Information

| Size | Cat. # | Shelf Life | Storage Conditions |
|-------|----------|--------------------------------------|------------------------|
| 500mL | APBS-500 | 12 months from date of manufacturing | Store dark at +15-30°C |